

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-072-0440	Revision: 13
Remarks		single use only

Specification

Product	Reservoir PS, 50ml, high profile, 1 bottom, universal purpose, white, grad.
Dimensions	length: 146,1 mm (+/- 0,2 mm)
	width: 76,2 mm (+/- 0,2 mm)
	height: 27,4 mm (+/- 0,2 mm)
Colour	white
Graduation	grad.
Nominal Capacity	50ml
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months
Packaging	pcs per bag: 5
	pcs per carton: 100
	pcs per pallet: 5.500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

