

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	06-021-0251	Revision: 2
Remarks	Suitable for various spectrophotometer. Each box selected from same mold cavity.	single use only

Specification

Product	Macro cuvette PS, 4ml, for photometer, cavity sorted, application 340-800 nm		
Dimensions	outer diameter: 12,50 mm (+/- 0,2 mm)		
	inner diameter: 10,00 mm (+/- 0,2 mm)		
	height: 45,00 mm (+/- 0,2 mm)		
Nominal Capacity	4ml		
Material Characteristics	polystyrene (PS)		
Temperature range	-50°C - 55°C		
Product Characteristics	biol. inert		
	physiol. harmless		
	no heavy metals acc. to EC/94/62		
	Food contact: raw material in compl. with reg. EC 1935/2004		
Shelf - Life	60 months at controlled storage conditions		
Packaging	pcs per carton: 1.000		
	pcs per pallet: 54.000		
	pcs per dsp case: 100		
Storage	dry condition: 40%-60% rel. humidity		
	bags / cases orig. closed		
	storage temp. +10°C to max. +30°C		
	pal. not stackable		
Approved by nerbe plus quality department			