

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-401-5003	Revision: 6
Remarks		single use only

Specification

Product	Inoculation loop PS, 1µl, flexible, 20 pcs Zip Bag green, sterile R
Dimensions	diameter loop: 0,90 mm (+/- 0,1 mm)
	diameter handle: 4,00 mm (+/- 0,1 mm)
	length: 199,00 mm (+/- 0,2 mm)
Colour	green
Nominal Capacity	1µl
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months
Packaging	pcs per inner carton: 1.000
	pcs per carton: 5.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	