

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 09-516-5009 | Revision: 10 |
| Remarks | | single use only |

Specification

| | | | | | | | | | | | | | |
|---|--|---------------------------|--------------------------------|---------------------------------------|---------------------------------------|-------|-------------------------------|---------------------------------------|--------|--------------------------------|--------------------------------|------------------------------------|---|
| Product | Transport swab with viscose tip, with aluminum stick and viscose tip, sterile R, Class 1s / MD CE2797 | | | | | | | | | | | | |
| Classification | Class 1s / MD CE2797 | | | | | | | | | | | | |
| Dimensions | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>tube:</td></tr> <tr><td>length: 145,00 mm (+/- 0,2 mm)</td></tr> <tr><td>outer diameter: 12,00 mm (+/- 0,2 mm)</td></tr> <tr><td>inner diameter: 10,70 mm (+/- 0,2 mm)</td></tr> <tr><td>plug:</td></tr> <tr><td>length: 20,00 mm (+/- 0,2 mm)</td></tr> <tr><td>outer diameter: 12,00 mm (+/- 0,2 mm)</td></tr> <tr><td>stick:</td></tr> <tr><td>length: 150,00 mm (+/- 0,2 mm)</td></tr> <tr><td>diameter: 0,90 mm (+/- 0,2 mm)</td></tr> <tr><td>diameter tip: 1,50 mm (+/- 0,2 mm)</td></tr> <tr><td>total length of tube and plug: 165,00 mm (+/- 0,2 mm)</td></tr> </table> | tube: | length: 145,00 mm (+/- 0,2 mm) | outer diameter: 12,00 mm (+/- 0,2 mm) | inner diameter: 10,70 mm (+/- 0,2 mm) | plug: | length: 20,00 mm (+/- 0,2 mm) | outer diameter: 12,00 mm (+/- 0,2 mm) | stick: | length: 150,00 mm (+/- 0,2 mm) | diameter: 0,90 mm (+/- 0,2 mm) | diameter tip: 1,50 mm (+/- 0,2 mm) | total length of tube and plug: 165,00 mm (+/- 0,2 mm) |
| tube: | | | | | | | | | | | | | |
| length: 145,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| outer diameter: 12,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| inner diameter: 10,70 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| plug: | | | | | | | | | | | | | |
| length: 20,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| outer diameter: 12,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| stick: | | | | | | | | | | | | | |
| length: 150,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| diameter: 0,90 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| diameter tip: 1,50 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| total length of tube and plug: 165,00 mm (+/- 0,2 mm) | | | | | | | | | | | | | |
| Material Characteristics | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>tube: polypropylene (PP)</td></tr> <tr><td>cap: polyethylene (PE)</td></tr> <tr><td>swab: aluminium (AL)</td></tr> <tr><td>tip of swab: viscose (VI)</td></tr> </table> | tube: polypropylene (PP) | cap: polyethylene (PE) | swab: aluminium (AL) | tip of swab: viscose (VI) | | | | | | | | |
| tube: polypropylene (PP) | | | | | | | | | | | | | |
| cap: polyethylene (PE) | | | | | | | | | | | | | |
| swab: aluminium (AL) | | | | | | | | | | | | | |
| tip of swab: viscose (VI) | | | | | | | | | | | | | |
| Product Characteristics | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>biologically inert</td></tr> <tr><td>physiologically harmless</td></tr> </table> | biologically inert | physiologically harmless | | | | | | | | | | |
| biologically inert | | | | | | | | | | | | | |
| physiologically harmless | | | | | | | | | | | | | |
| Purity Grade | manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1) | | | | | | | | | | | | |
| Sterilisation Process | radiation sterilization SAL 10-6 | | | | | | | | | | | | |
| Shelf - Life | 48 months | | | | | | | | | | | | |
| Packaging | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>pcs per inner carton: 500</td></tr> <tr><td>pcs per carton: 2.000</td></tr> <tr><td>pcs per pallet: 48.000</td></tr> </table> | pcs per inner carton: 500 | pcs per carton: 2.000 | pcs per pallet: 48.000 | | | | | | | | | |
| pcs per inner carton: 500 | | | | | | | | | | | | | |
| pcs per carton: 2.000 | | | | | | | | | | | | | |
| pcs per pallet: 48.000 | | | | | | | | | | | | | |



Technical Datasheet

| | |
|--|--------------------------------------|
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +2°C to max. +30°C |
| | pal. not stackable |
| Approved by nerbe plus quality department | |