

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-023-0100	Revision: 1
Remarks		single use only

Specification

Product	Screw cap HDPE, for mailing container 30x85 mm and 31x126 mm, natural
Colour	natural
Material Characteristics	high density polyethylene (HDPE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 54.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	