

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-162-0100	Revision: 1
Remarks		single use only

Specification

Product	Mailing container PP, 60ml, Ø30x126 mm, with absorbing inlay, screw thread and print according to EN 829, without cap, natural, autoclavable up to 121°C
Dimensions	diameter: 30,00 mm (+/- 0,2 mm)
	height: 126,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	60ml
Material Characteristics	polypropylene (PP)
Temperature range	-196°C - 121°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 300
	pcs per carton: 300
	pcs per pallet: 7.200
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

