

# Technical Datasheet

Date: 13.10.23

Page 1 of 1

|                        |             |                 |
|------------------------|-------------|-----------------|
| EN ISO 13485 certified |             |                 |
| Article - No:          | 11-222-0100 | Revision: 1     |
| Remarks                |             | single use only |

## Specification

|                                           |                                                                                                                                                          |
|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product                                   | Mailing container PP, 50ml, Ø44x114 mm, with absorbing inlay, screw thread and print according to EN 829, without cap, natural, autoclavable up to 121°C |
| Dimensions                                | diameter: 44,00 mm (+/- 0,2 mm)                                                                                                                          |
|                                           | height: 114,00 mm (+/- 0,2 mm)                                                                                                                           |
| Colour                                    | natural                                                                                                                                                  |
| Nominal Capacity                          | 50ml                                                                                                                                                     |
| Material Characteristics                  | container: polypropylene (PP)                                                                                                                            |
|                                           | inlay: cellulose                                                                                                                                         |
| Temperature range                         | -196°C - 121°C                                                                                                                                           |
| Autoclavability                           | autocl. up to 121°C                                                                                                                                      |
| Product Characteristics                   | biologically inert                                                                                                                                       |
|                                           | physiologically harmless                                                                                                                                 |
|                                           | no heavy metals according to EC/94/62                                                                                                                    |
| Packaging                                 | pcs per bag: 250                                                                                                                                         |
|                                           | pcs per carton: 250                                                                                                                                      |
|                                           | pcs per pallet: 5.500                                                                                                                                    |
| Storage                                   | dry condition: 40%-60% rel. humidity                                                                                                                     |
|                                           | bags / cases orig. closed                                                                                                                                |
|                                           | storage temp. +10°C to max. +30°C                                                                                                                        |
|                                           | pal. not stackable                                                                                                                                       |
| Approved by nerbe plus quality department |                                                                                                                                                          |