

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-802-6100	Revision: 1
Remarks		single use only

Specification

Product	Multi purpose container PP, 200ml, Ø66x88 mm, with fitted snap-on lid PE, transparent, graduated, sterile R
Dimensions	outer diameter: 66,00 mm (+/- 0,2 mm)
	inner diameter: 62,00 mm (+/- 0,2 mm)
	height (with lid): 88,00 mm (+/- 0,2 mm)
	height (without lid): 87,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	200ml
Temperature range	-80°C - 75°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 220
	pcs per carton: 220
	pcs per pallet: 4.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department