

Technical Datasheet

Date: 13.10.23

Page 1 of 1

EN ISO 13485 certified		
Article - No:	11-853-0100	Revision: 1
Remarks		single use only

Specification

Product	Lid PE for multi purpose container 11-852-0100, white
Colour	white
Material Characteristics	polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 3.750
	pcs per carton: 3.750
	pcs per pallet: 540.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	