

Technical Datasheet

Date: 18.05.26

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-001-5000	Revision: 6
Remarks	Universal fit	single use only

Specification

Product	GTS-pipette, Single-channel, variable, 500-5.000µl color code purple
Colour	color code purple
Nominal Capacity	500-5.000µl
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Shelf - Life	3 years warranty
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	