

# Technical Datasheet

Date: 18.05.26

Page 1 of 1

|                        |               |                 |
|------------------------|---------------|-----------------|
| EN ISO 13485 certified |               |                 |
| Article - No:          | 12-008-0100   | Revision: 6     |
| Remarks                | Universal fit | single use only |

## Specification

|                                           |                                                                           |
|-------------------------------------------|---------------------------------------------------------------------------|
| Product                                   | GTS-pipette, 8-channel, variable, 10-100µl, color code yellow             |
| Colour                                    | color code yellow                                                         |
| Nominal Capacity                          | 10-100µl                                                                  |
| Purity Grade                              | manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1) |
| Shelf - Life                              | 3 years warranty                                                          |
| Storage                                   | dry condition: 40%-60% rel. humidity                                      |
|                                           | bags / cases orig. closed                                                 |
|                                           | storage temp. +10°C to max. +30°C                                         |
|                                           | pal. not stackable                                                        |
| Approved by nerbe plus quality department |                                                                           |