

Technical Datasheet

Date: 06.11.25

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-012-0010	Revision: 6
Remarks	Universal fit	single use only

Specification

Product	np GTS-pipette, 12-channel, variable, 0,5-10µl, color code grey
Colour	color code grey
Nominal Capacity	0,5-10µl
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	