

### Electrodes for ORP measurement

2000233	Standard ORP	electrode with platinum ring for ORP/Redox
2000663	Plast ORP	electrode with platinum pin for ORP/Redox
2000673	Plast ORP/BNC	electrode with platinum pin for ORP/Redox, fix cable BNC

### Conductivity cells

3004062	VP8	conductivity cell, cell constant c=1, fix cable BNC
3003352	VPT 80/1	conductivity cell, cell constant c=1, fix cable BNC, temperature sensor
3004012	VPT51/01	conductivity cell, cell constant c=0,1, fix cable BNC, temperature sensor for ultrapure water
3004072	VPT 80/10	conductivity cell, cell constant c=10, fix cable BNC, temperature sensor for sea water
3004002	2301T	conductivity cell, cell constant c=1, fix cable BNC, temperature sensor

### Buffer solution

800033	Buffer pH4, $\pm 0,01$ / 25°C, 250ml bottle, red
800133	Buffer pH4, $\pm 0,01$ / 25°C, 500ml bottle, red
800233	Buffer pH4, $\pm 0,01$ / 25°C, 5 Liter red
800043	Buffer pH7, $\pm 0,01$ / 25°C, 250ml bottle, green
800143	Buffer pH7, $\pm 0,01$ / 25°C, 500ml bottle, green
800243	Buffer pH7, $\pm 0,01$ / 25°C, 5 Liter green
800073	Buffer pH9,21, $\pm 0,02$ / 25°C, 250ml bottle, blue
800173	Buffer pH9,21, $\pm 0,02$ / 25°C, 500ml bottle, blue
800273	Buffer pH9,21, $\pm 0,02$ / 25°C, 5 Liter blue
800063	Buffer pH10, $\pm 0,02$ / 25°C, 250ml bottle, colorless
800163	Buffer pH10, $\pm 0,02$ / 25°C, 500ml bottle, colorless
800263	Buffer pH10, $\pm 0,02$ / 25°C, 5 Liter colorless



Buffer solutions

*DAKKS certified buffer solutions are available as well*

800613	Conductivity standard 84 $\mu\text{S}/\text{cm} \pm 1\%$ / 25°C, 500 ml bottle
800633	Conductivity standard 1.413 $\mu\text{S}/\text{cm} \pm 1\%$ / 25°C, 500 ml bottle
800643	Conductivity standard 12.880 $\mu\text{S}/\text{cm} \pm 1\%$ / 25°C, 500 ml bottle
850513	Conductivity-Standard 5 $\mu\text{S}/\text{cm} \pm 0,1$ / 25°C, 300 ml bottle with N.I.S.T Certificate
850313	Conductivity-Standard 5 $\mu\text{S}/\text{cm} \pm 0,1$ / 25°C, 300 ml bottle with DFM Certificate
850303	Conductivity-Standard 1,3 $\mu\text{S}/\text{cm} \pm 0,1$ / 25°C, 300 ml bottle with DFM Certificate

### Redox- buffer

800303	250 ml Redox/ORP 200 mV $\pm 5\text{mV}$ / 25°C
800313	250 ml Redox/ORP 475 mV $\pm 5\text{mV}$ / 25°C
800323	500 ml Redox/ORP 475 mV $\pm 5\text{mV}$ / 25°C



Buffer solutions



Cond-Solutions



Printer 1001502

### More accessories

1001733	Cable S7 / BNC 1 meter
1001743	Cable S7 / DIN 1 meter
1002012	Temperature sensor NT55
1000112	Elektrode holder
1001502	Printer for EC-35 and EC-45